

AXIAL FLOW PUMP TEST RIG MAKE – LABTEK, ORIGIN – INDIA



Pump : Axial Flow Pump to Discharge About 150LPM.

Motor : 3HP, Three phase. 400/440V, AC supply.

Electrical Panel : Switch & Energy meter for the measurement of power Input.

Sump Tank : M.S. Unit with FRP Lining/S.S Unit of suitable size to store sufficient water for independent circulation through the unit.

Measuring Tank : M.S. Unit with FRP Lining/ S.S. Unit of suitable size provided with standard fittings for measuring discharge of water.

Others : Piping system consisting of Pipes, Pressure relief valve, Gate valve and fitting with pressure gauge.



Sun LabTek Equipments (I) Pvt. Ltd.
ISO:9001:2008 Certified Company

