

Digital Constant Water Bath

MAKE – LABTEK, ORIGIN – INDIA



Same as ZI 11010 but supplied with Digital Temperature Controller range ambient to 100.0 deg C accuracy $\pm 0.2^{\circ}$ C, least count $\pm 0.1^{\circ}$ C. & FHP Motor cum Stirrer. Heaters are provided according to the size of the chamber. A control panel in which the digital controller, switches for mains, FHP stirrer, regulator for the stirrer is fit is by the side. Sizes available: 18"x12"x10" inner chamber in stainless steel, outer in mild steel. Rating of heater 2000 W suitable to run on 230 V A.C. Single Phase.