

# Digital Double Door Drying Oven

MAKE – LABTEK, ORIGIN – INDIA



Digital Double Door Drying Oven

This range of ovens is designed for drying large quantities of soil and aggregate samples and maintains temperature in accordance with most B.S. and ASTM requirements. The 225 Lts. and 425 Lts.

## OPTIONS AND ACCESSORIES :

Dial thermometer with collar fixing for door Heavy duty shelves (on Request at an extra price) Color of the Oven (Only on qty. greater than 3 Units)

## Specification :

1. The larger soil drying ovens are constructed with mild steel and powder coated exterior and Stainless Steel coated chamber which is both durable and corrosion resistant.
2. All Sizes of Oven are fitted with a (PID Controller) having a microprocessor digital control as standard.
3. All units comes with MCB tripping device in cases of any short circuit, overloading of voltage and Low voltage cut off. Wiring of heaters can be connected / dis-connected removing the side panels of the Oven.