

Humidity Chamber

MAKE – LABTEK, ORIGIN – INDIA



CONSTRUCTION :

Double walled outer M.S. Sheet duly powder coated. Inner Aluminium/ Stainless Steel with 2nos. adjustable shelves duly perforated. Having heavy duty M.S. Stand. The door is an impressive structure having glass window & sealed by Rubber Gasket. Gap between the walls is filled with glass.

CONTROL PANEL :

Consist of Temp. controller Le. Thermostat for dry temp. Control, Wet Temp. controller Energy Regulator, Mains, Dry, Wet Indicator Humidity Indicator, Mains, Illumination Switch. Humidity controlled by Imported Humidistate placed inside the chamber. Supplied with Hygrometer/ wet and dry thermometer to check exact Humidity inside the chamber.