

Laboratory Incubator

MAKE – LABTEK, ORIGIN – INDIA



Incubators are designed fabricated & tested to suit various applications in growing field of Medical Agriculture. Research Laboratories of Hospital/Bio-Tech lab. Suitable to work on 230 V A.C. Supply 50 Hz

Construction :

Double walled construction outer body made of M.S. Duly stove enamel/powder coated, with inside Aluminium/Stainless Steel Chamber having 60-75mm glass wool insulation, door with glass window

HEATING :

Heating element are made from high grade (Kanthal) Resistant wire. Heaters are fitted/wound along all sides & bottom. Complete with 2 pilot lamp cord & plug pin etc.

OPTIONAL :

1. Air circulation can be fitted in all the above size at extra xost.
2. Digital Temp. Controller/Indicator can be fitted with air circulation system at system at extra cost if required.