

## LOVIBOND TINTOMETER MAKE – LABTEK, ORIGIN – INDIA



The Comparator System 2000 is based on a versatile 2-field visual comparator with an optional lighting unit to guarantee correct lighting conditions for colour grading and colorimetric analysis.

Version	Application	Scales
Model F	Edible Oils & Fats, Food & Beverages, Petroleum Fuels, Oils & Waxes, Pharmaceuticals & Cosmetics, Beers, Malts & Caramels	Lovibond® RYBN Colour, Kreis Value (Rancidity in Fats)
Model F (BS 684), with compensating slides	Edible Oils & Fats, Food & Beverages, Petroleum Fuels, Oils & Waxes	BS 684