

Laboratory Incubator

MAKE – LABTEK, ORIGIN – INDIA



Seed Germinator

Most suitable for conducting various experiments on a variety of seeds under different conditions of temperature and humidity. Atmospheric conditions are created within the chamber to perform different chamber made up of anodized Aluminium or S.S. and outer chamber of M.S. Sheet duly painted. Angle Iron frame with caster wheel at bottom for easy movement. Humidity (not adjustable) is created in a S.S. water reservoir fitted at the bottom to provided 90% to 95% $\hat{\pm}3\%$ RH Humidity.

OPTIONAL ACCESSORIES :

1. Exterior lighting arrangement (on both sides)
2. Automatic Timer (0-24 hours)
3. Electronic Digital Temperature Indicator-cum-Controller.
4. Dual Setting Digital Temperature controller for setting two different Temperature for day and night.

Inside Chamber Size (WxHxD) : 555mm x 910mm X 605mm LTR.: 308