

TILTING FLUME APPARATUS MAKE – LABTEK, ORIGIN – INDIA



TECHNICAL SPECIFICATIONS:

A 76mm wide, 250mm high open channel.

Clear acrylic sides to give visibility of the working section.

Ajacking system permits the slope of the channel bed to be adjusted between -1% and +3%.

Service module: Made from Acrylic

Sump tank capacity: 250 litres (S. S.)

Volumetric tank capacity: 40 litres (GRP)

Inlet tank with flow stilling arrangement.(GRP)