

Universal Water Bath

MAKE – LABTEK, ORIGIN – INDIA



Double walled inner chamber in stainless steel outer chamber in mild steel temperature range Ambient to 100° C controlled by capillary thermostat. Accuracy $\pm 1^{\circ}$ C operated on 230V A.C. Single Phase. These baths are suitable for washing Marshall specimen. Lid is made up of stainless steel with a hole for glass thermometer to penetrate.

Following sizes are available:

- A) 12"lengthx 9" width x 7" height having rating 750 W.
- B) 18"length x 12" width x 7" height having rating 1500 W.

Customised sizes also available.