

AIRBORNE PARTICLE COUNTER

MODEL NO SL -EE -019, MAKE - LABTEK, ORIGIN - INDIA



Functions:

1. Eight particle-size channels (0.3, 0.5, 0.7, 1.0, 2.0, 5.0, 7.0, 10 μ m) 9-digit LED display (7 digitals for data display)
2. Count data is displayed as cumulative or differential count for each particle-size channel
3. Temperature and humidity display (optional)
4. RS-232 communication interface
5. Date and clock display
6. 95%UCL calculation
7. Built-in printer, multi-print mode
8. Data storage 200 samples