

ANGLE DEKKOR

MODEL NO – SL-MEL-19, MAKE – LABTEK, ORIGIN – INDIA



- a). Micro optic and manual operated.
- b). Maximum working distance : 12m
- c). Accuracy : 6 sec. of arc
- d). Resolution: 0.5 sec. of arc

FUNCTIONS:

- a). Measures surface flatness using collimated light and a moving lens to focus the beam as it travels across the test surface.
- b). Any angle can be measured by comparison with a standard angle or any change of angle can be observed instantly and accurately.