

## COMPACTION PEDESTAL/COMPACTION HAMMER

MODEL NO – SL-GT-018/19, MAKE – LABTEK, ORIGIN – INDIA



### Compaction Pedestal

Comprising a 12" square x 1 inch thick steel plate secured to an 8" square x 18" high wooden block. 4 angle brackets are supplied for securing the block. A specimen mould holder is fitted to the steel plate.

### Compaction Hammer

Satisfies BS 598. The hammer has a 4535 g sliding weight with a free fall of 457mm.