

PERMEABILITY APP. (3 CELL MODEL) MODEL NO SL -CC -101, MAKE - LABTEK, ORIGIN - INDIA



As Per IS 3085 / DIN Standard (Single Cell Model)

One of the durability test of concrete is to determine permeability of water through specimen. Permeability apparatus is used for determining the permeability of cement mortar and concrete specimens of 15cm cubes cast in the laboratory.

SPECIFICATION :

The concrete permeability apparatus comprises of a brass / gunmetal cell of Square/Round cross-section mounted on a stand and a pressure chamber is connected to the cell through copper tubing and T-connector mounted on the stand with sleeve packed valve and rubber hose pipe with end connections.