

CONDUCTIVITY METER

MODEL NO SL –EE –009,010,011, MAKE – LABTEK, ORIGIN – INDIA



DIGITAL PH METER (MODEL SL-EE-006)

MICROPROCESSOR CONDUCTIVITY / TDS / TEMP. (MODEL SL-EE-009)

DISPLAY	16 X 2 Alpha Numeric LCD Display
RANGE	0.1 μ S to 100mS (6 Decadic Ranges)
RESOLUTION	0.01 μ S.
ACCURACY	$\hat{A}\pm 1\%$ of F.S. $\hat{A}\pm 1$ Digit
CALIBRATION	Auto / Manual.
TEMPERATURE COMP.	0 to 100 UC.
CELL CONSTANT	0.1 To 5.0 Adjustable.
SENSOR	Up to 100 Samples.
POWER	230V $\hat{A}\pm 50$ Hz. AC.

DELUXE DIGITAL CONDUCTIVITY METER (MODEL SL-EE-010)