

## CONDUCTIVITY-TDS METER

MODEL NO SL -EE -012, MAKE - LABTEK, ORIGIN - INDIA



### SPECIFICATIONS

#### SPECIFICATIONS

<b>Range</b>	
Conductivity	2 $\mu$ S To 200mS With 0.1 $\mu$ Cell Constant And 20 $\mu$ S To 200mS With Unity Cell Constant in 5 Decade Range
TDS	2 PPM To 2 PPT With 0.1 Cell Constant And 20 PPM To 200 PPT With Unity Cell Constant, in 5 Decade Range
Temperature	0 $^{\circ}$ C To 100 $^{\circ}$ C

#### Resolution