

CONSTANT VOLUME MOULD

MODEL NO SL -GT - 036, MAKE - LABTEK, ORIGIN - INDIA



Often in the laboratories it is required to quickly prepare remoulded soil specimen for Triaxial and Unconfined tests. Constants volume moulds are used for this purpose, using either dynamic or static compaction. The moulds are available in different sizes. For specimen size 38 mm dia x 76 mm long comprises. Split mould 38 mm dia x 126 mm long. End plugs 38 mm dia x 25 mm long pair. Split collar interchangeable. Ejecting plunger 38 mm dia x 126 mm long.

Accessories split mould 38 mm dia x 76 mm long and Split mould 38 mm dia x 86 mm long.