


DIGITAL DISSOLVED OXYGEN MODEL NO SL -EE -016,017, MAKE - LABTEK, ORIGIN - INDIA



DIGITAL DISSOLVED OXYGEN ANALYZER CUM TEMP. METER (MODEL SL-EE-016)

RANGE	0 to 20.0ppm
RESOLUTION	0.1 ppm
ACCURACY	$\hat{A}\pm 0.2\text{ppm}$
TEMPERATURE COMP.	0to 50 deg. C. Manual.
DISPLAY	3.5 Digit LED
STIRRER	Provided.
POWER	230V $\hat{A}\pm 10\%$ 50 Hz
SENSOR	DO Amperometric Gold / Sliver Temp. RTD (PT-100)
TIME RESPONSE	10 Sec

DIGITAL DISSOLVED OXYGEN METER (MODEL SL-EE-017)