

DIGITAL PH METER

MODEL NO SL –EE –006,007,008, MAKE – LABTEK, ORIGIN – INDIA



DIGITAL PH METER (MODEL SL-EE-006)

DIGITAL PH METER (MODEL SL-EE-006)

SPECIFICATIONS

PH Range	0 To 14 PH
MV Range	0 To +1999 MV
Resolution	0.01 PH, In PH Range 1 Mv In MV Range
Repeatability	0.01 PH 1 Digit 1 MV 1 Digit
Standardization Range	Min. 2 PH.
Temp Compensation	0 - 1000c
Display	3½ Digits LED
Power	230v 10%, 50 Hz

DIGITAL PH METER (MODEL SL-EE-007)

B-3A, Shiv Shakti Complex, East Vinod Nagar, 1st Floor, New Delhi - 110 091

Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915

E-mail: info@labtekindia.com, Web: www.labtekindia.com