


# SOFTWARE FOR TRIAXIAL TEST (ALL ELECTRONIC MODELS)

MAKE – LABTEK, ORIGIN – INDIA



- 1). Does all calculations of UC, UU, CU, CU & CD triaxial tests
- 2). Has option for manual as well as automatic recording of data
- 3). Displays the following plots (Graphical)
  - a). Plots consolidation curve & evaluates  $t_{100}$  and calculates strain rate depending Upon the drainage condition in CD & CU
  - b). Stress - strain curves after every test
  - c). Plots pore water pressure and A - factor (for CU) and volumetric strain (for CD) vs. axia Effective stress ratio vs. axial strain
  - d). p-q plot
  - e). Mohr stress circles & envelope giving C & in terms of effective stress.