

## FRANCIS TURBINE TEST RIG

MAKE – LABTEK, ORIGIN – INDIA



## FRANCIS TURBINE TEST RIG

### Francis Turbine

<b>Type</b>	: Inward Flow Reaction Turbine
<b>Capacity</b>	: 5HP
<b>Rated Speed</b>	: 1250 RPM
<b>Discharge Capacity</b>	: 1800 Liters/Minute
<b>Supply Head</b>	: 25 Meters
<b>Guide Vanes</b>	: Gun Metal Vanes (Aerofoil Blade Shaped)
<b>Loading</b>	: Rope Brake
<b>Material</b>	: Cast Iron
<b>Drum Size</b>	: 300 mm diameter

### Supply Pump Set