

HARDNESS TESTER FOR MASTIC ASPHALT MODEL NO SL -BT -006, MAKE - LABTEK, ORIGIN - INDIA



IS 1195

For determining the hardness number of Mastic Asphalt for flooring. It consists of an internally insulated cabinet to the base of which is fitted a water bath having two taps. The bath is heated by an immersion heater and the temperature is controlled with the thermostat at $35\text{Å}^{\circ} \pm 0.5\text{Å}^{\circ} \text{C}$. A 6.5mm dia pin is loaded on the specimen by means of a lever to give a 31.7 Kg weight. A dial gauge 0.01 x 25 mm is provided to record the penetration of pin into the specimen.