

HEAT TRANSFER THROUGH LAGGED PIPE

MAKE – LABTEK, ORIGIN – INDIA



01. Pipes

a. Diameter	: 50mm/150mm Inner/Outer (Approx.)
b. Length	: 450mm (Approx.)
c. Material	: G.I.
d Quantity	: 1 No. Each

02. Measuring Instruments

a. Digital Voltmeter Range	: 0 to 300 V
b. Digital Ammeter & Range	: 0 to 5 Amps

c. Digital Temperature Indicator with Selector switch