

## MORTAR NEEDLE PENETROMETER

MODEL NO – SL-CC-098, MAKE – LABTEK, ORIGIN – INDIA



It is used for finding out the rate of hardening of mortar sieved from concrete spring and a stem graduated from 0-70 kg x 1 kg. Six interchangeable penetration needles of areas 645, 323, 65 32 and 16mm sq. Is provided . The penetration resistance is measured by the force exerted to penetrate the mortar by 25mm and is indicated by a sliding ring on the stem, which is graduated. Needle shanks are marked at every 12.5mm. Complete in a wooden carrying case.

**B-3A, Shiv Shakti Complex, East Vinod Nagar, 1<sup>st</sup> Floor, New Delhi - 110 091**

**Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915**

**E-mail: info@labtekindia.com, Web: www.labtekindia.com**