

## RIFFLE SAMPLE DIVIDER

MODEL NO SL -BT -006, MAKE - LABTEK, ORIGIN - INDIA



**As Per IS: 1607-1960**

**Specification :**

Riffle sample divider is required for the rapid collection of true representative samples from aggregates, sand and fillers. Riffle Sample Divider consists of a metal box fitted with a series of chutes of equal width which discharge the material alternatively in opposite directions into separate pans. The chutes of the riffle are steep enough to allow rapid flow of the material. The sample divider is supplied with two containers and a sampling scoop.