

## 'ROTAP' SIEVE SHAKER

MODEL NO SL -GT - 088, 089, MAKE - LABTEK, ORIGIN - INDIA



### 'ROTAP' Sieve Shaker (SL-GT-088)

To make process of sieving simpler and quicker, Rotap sieve shaker is useful , this produces circular shaking of the sieves. At the same time the sieves are tapped. The mechanism for imparting circular action and topping is oil immersed in a tank and in motorized. The shaker can accommodate upto 7 nos. Sieves of dia 200 or 300mm, 450mm Suitable operation of 230 V A.C. Single Phase.