

SEMI-AUTOMATIC PENETROMETER MODEL NO SL -BT -005, MAKE - LABTEK, ORIGIN - INDIA



IS 310, 1203, 1448, IP 60, 49, 50 ASTM D5, 217, D637, BS 2000-49

Same as above but the unit is compact with timer to control duration of penetration. The instrument is provided with lead screw gear arrangement, Leveling screws, Spirit level.

B-3A, Shiv Shakti Complex, East Vinod Nagar, 1st Floor, New Delhi - 110 091

Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915

E-mail: info@labtekindia.com, Web: www.labtekindia.com