

SPRING TESTING MACHINE MAKE – LABTEK, ORIGIN – INDIA



TECHNICAL SPECIFICATION

- Maximum Capacity 500 kg / 5000N
- Type of Spring Tension & Compression
- Free length of tension spring 300mm
- Diameter of compression spring 200mm
- Load Application By top hand wheel
- Maximum Cross Head Travel-without sample/grip 500mm
 - Force Measuring Resolution 0.1kg(100gms) / 1N
 - Pitch 5mm

One complete revolution of Handle (360 degree) Travel of Crosshead up or down

- Type of Construction Double Column
- Elongation / Deflection / Least Count 1 ft scale(300mm)

B-3A, Shiv Shakti Complex, East Vinod Nagar, 1st Floor, New Delhi - 110 091

Phone : 011-65708108, 65808108 Mobile No : 09971077233, 09811888915

E-mail: info@labtekindia.com, Web: www.labtekindia.com