

## Universal Testing Machine

MAKE – LABTEK, ORIGIN – INDIA



**Loading Accuracy:** As high as + 1 %

**Speeds:** Straining at variable speeds to suit wide range of materials

**Facilities for Tests:** Motor-driven threaded columns for quick and convenient adjustment of lower crosshead to facilitate rapid fixing of test specimen.

**Autographic Recorder:** Simultaneous roll autographic recorder supplied as standard to enable study of the behaviors of materials.

**Ideal Dia:** High reading accuracy due to large size ideal design of dial.

**Large Columns:** Large effective clearance between columns enables testing of standard specimen as well as structures



**Sun LabTek Equipments (I) Pvt. Ltd.**  
ISO:9001:2008 Certified Company

