

## MORTAR CUBE VIBRATING MACHINE

MODEL NO – SL-CC-092, MAKE – LABTEK, ORIGIN – INDIA



Vibrating Machine is used for vibrating the mix in moulds at a frequency of  $12,000 \pm 400$  cycles per minute. The vibrator is mounted over 4 coiled springs and the vibrations are developed by means of a revolving eccentric shaft. The centre of gravity of the vibrator, including the cube mould, is either at the centre of eccentric shaft or within 25mm below it. The simple design of the machine facilitates easy assembly and dismantling of the cube moulds.