

## WET SIEVE SHAKER, HAND OPERATED MODEL NO SL -GT - 090, 091, MAKE - LABTEK, ORIGIN - INDIA



### **Wet Sieve Shaker, Hand Operated (SL-GT-090)**

It consists of a motor with reduction gear box mounted at the top of the reservoir. The sieve holder, capable of accommodating seven sieves of 150 mm or 200 mm diameter is hung from an eccentric which imparts up and down motion to the nest of sieves held in the sieve holder and immersed in water. Supplied without sieves, and lid & receiver.