

## AGGREGATE IMPACT VALUE TESTER

MODEL NO SL -BT -008, MAKE - LABTEK, ORIGIN - INDIA



### As Per IS 2386 (P-IV)

#### Specification :

Aggregate Impact Value Test Apparatus is used for determining the aggregate impact value. The Sturdy Construction consists of a base and support columns to form a rigid framework around the quick release trigger mechanism to ensure an effective free fall of the hammer during test. The free fall can be adjusted through  $380\pm 5\text{mm}$ . The hammer is provided with a locking arrangement. Aggregate impact value test apparatus is supplied complete with a cylindrical measure of 75mm dia x 50mm depth, an automatic blow counter and a tamping rod.