

DIGITAL TURBIDITY METER

MODEL NO SL -EE -014,015, MAKE - LABTEK, ORIGIN - INDIA



SPECIFICATIONS

DIGITAL TURBIDITY METER (MODEL SL-EE-014)

DISPLAY	3.5 Digital LED.
POWER	230v $\hat{\pm}$ 10% 50 Hz.
RANGE	0 to 200 NTU.
RESOLUTION	1 NTU.
ACCURACY	$\hat{\pm}$ 2% FS $\hat{\pm}$ 1 Digit.
SAMPLE SYSTEM	30mm Clear Glass Test Tubes.
LIGHT SOURCE	6V-1Amp. Tungsten Lamp.
DETECTOR	Photocell / Photodiode.

DIGITAL TURBIDITY METER (MODEL SL-EE-015)