

## CUT SECTION MODEL OF SEMI FLOATING DEFERENTIAL AND WHEEL MECHANISM (WORKING) MODEL NO – SL -80012, MAKE – LABTEK, ORIGIN – INDIA



This model is made out of used BAJAJ engine, suitably sectioned to show the internal construction of the engine, gear box, clutch and rear wheel mechanism and the model is fitted on a sturdy iron frame. The complete internal details can be demonstrated by operating the kicker lever to show the working of the piston, sparks from the spark plugs can also be shown